## Le Ballon Verde, alcohol breathalyzer for single use - 0.2 ‰

This breathalyzer for single-use is calibrated at 0.2 promille Blood Alcohol Concentration (BAC). The breathalyzer is NF certified and is REACH compliant so complies to European legislation. The test can be delivered in standard Green packaging or with personalized label.

Le Ballon Verde is easy to use: remove the test tube from the edge of the balloon, inflate Le Ballon, compress the ends of the tube and push the air in the balloon through the tube. Within a minute you can read the result. If the white colored crystals turn red, you have exceeded the permitted limit of 0.2 per mille.

The alcohol test is supplied in standard French version. The alcohol test can be provided with instruction in other languages including Dutch, French, English and German in combination with a personalized label or as a single instruction per bulk packaging. The tester is calibrated to 0.2 g/L (permille) blood alcohol concentration.

### Characteristics

- One-time use
- NF X 20-702: Oct. 2014 certified
- THT of 24 months after delivery
- Hygienically packaged
- Compact and lightweight
- Reliable and accurate
- Complies with REACH 2017/09
- User-friendly
- Cheap, works without batteries and no calibration needed

### Specifications

|               |                             |
|---------------|-----------------------------|
| Dimensions    | 54x100x15mm                 |
| Weight        | 15gr                        |
| Certification | NF-X 20-702 of October 2014 |

|                           |  |
|---------------------------|--|
| Calibration               | 0.09 mg/l breath alcohol<br>0.20 g/L blood alcohol (promille, ‰)   |
| Environment               | 100% recyclable,<br>REACH compliant,<br>may be disposed as household waste                               |
| Result                    | Crystals remain white: your BAC is below 0.2 ‰<br>Crystals turn red: your BAC is equal to or over 0.2 ‰. |
| Accuracy                  | 100% accurate according to the test requirements of the NF standard                                      |
| Instruction               | French<br>Other languages available on request   |
| Shelf life                | 24 months  |
| Personalization possible: | Yes  |

## Manual

### Instructions For Single-Use ALCOHOL BREATHALYSER

Wait at least 10 minutes after drinking alcohol or other substances (mouthwash, medicine...) before you take the breath test.  
The maximum blood alcohol level is reached about 45 min to 60 min after last consuming alcohol.

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Push the ends in as far as you can.  
Enfoncer les embouts au maximum.

Blow up the bag with a single breath.  
Gonfler le ballon de façon continue et sans interruption.

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### Reading the result

: Open your hand and immediately check the colour of the crystals!

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